

What you need to do:

The detection efficiency calibration of real germanium detector arrays!

- Fit gamma-ray photopeaks with GRSISort
- Edit and run the RootEffi script
- Submit screenshots of your fitted peaks, edited RootEffi script and efficiency curve ©

What you need to know:

- How semiconductor detectors work
- How to find nuclear information
- How to fit peaks with GRSISort
- How to run a script in ROOT/GRSISort
- How to determine the detection efficiency





GRIFFIN, TRIUMF, Vancouver, Canada

64 crystals






MINIBALL, CERN, Geneva, Switzerland

24 crystals

GAMKA, iThemba LABS, Cape Town, South Africa

52 crystals





TIGRESS, TRIUMF, Vancouver, Canada


56 crystals



Configurations of Germanium Detectors

GRIFFIN, TRIUMF, Vancouver, Canada

64 crystals (16 clovers, not segmented)





Switzerland 24 crystals (8 triple clusters,

Closed-end coaxial right cylinders ©

GAMKA. iThemba LABS, Cape Town, South Africa

52 crystals (13 clovers, 4 segmented and 9 not segmented)






TIGRESS. TRIUMF, Vancouver, Canada

56 crystals (14 clovers, 8-fold segmented)


How to Determine the Detection Efficiency



Absolute Efficiency $\epsilon_{\gamma} = \frac{\text{Number of } \gamma\text{-rays detected}}{\text{Number of } \gamma\text{-rays emitted}}$

$$= \frac{N_{\gamma, detected}}{I_{\gamma} \cdot A \cdot t} \leftarrow$$

Area of the photopeak
Data collection time
Activity of the source
Absolute intensity of
the gamma ray



$$\epsilon_{\gamma}(E) = 10^{p_0 + p_1 \log(E) + p_2 \log^2(E) + p_3/E^2}$$

Relative Efficiency
$$\epsilon_{\gamma} = \frac{\frac{N_{\gamma 1, detected}}{I_{\gamma 1} \cdot A \cdot t}}{\frac{N_{\gamma 2, detected}}{I_{\gamma 2} \cdot A \cdot t}} = \frac{N_{\gamma 1, detected}/I_{\gamma 1}}{N_{\gamma 2, detected}/I_{\gamma 2}}$$

How to Find Nuclear Information

- Lund/LBNL Nuclear Data Search: http://nucleardata.nuclear.lu.se/toi/
- National Nuclear Data Center: https://www.nndc.bnl.gov/

The Lund/LBNL Nuclear Data Search

Version 2.0. February 1999

S.Y.F. Chu¹, L.P. Ekström^{1,2} and R.B. Firestone¹

¹ LBNL, Berkeley, USA
² Department of Physics, Lund University, Sweden


WWW Table of Radioactive Isotopes

Radiation search
Nuclide search
Atomic data (X-rays and Auger electrons, very preliminary!)
Periodic chart interface to the nuclides
Summary drawings for A=1-277 (PDF)
Nuclear charts (PDF, 333 kbyte)


Database status

Table of Isotopes (Tol)

About this service Tol home page




Semiconductors vs. Scintillators




Material Property	Nal(TI)	Germanium	Silicon
Туре	Scintillator	Semiconductor	Semiconductor
Cooling	No	~ 95 K	No
Density (g/cm²)	3.7	5.3	2.3
Band gap energy (eV)	20	2.9	3.8
Energy resolution at 1332 keV (keV)	13	1.9	

How Semiconductor Detectors Work

- An external field creates an area depleted of free charge carriers.
- Radiation interacts with the crystal and produces of electron-hole pairs.
- The electrons and holes drift towards electrodes and the electric pulse is amplified and processed.





K.S. Krane, *Introductory Nuclear Physics* (1988). G.F. Knoll, *Radiation Detection and Measurement* (1989).

W.R. Leo, Techniques for Nuclear and Particle Physics Experiments (1994).